

Human HSPH1 protein (Recombinant) (His-Tag) (STJP019665)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 50mM NaCl
Molecular Weight	100.9 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10808
Gene Symbol	HSPH1
UniProt ID	HS105_HUMAN
Immunogen Region	1-858aa
Immunogen Sequence	MSVV GLDVGSQLSCY IAVARAGGIE TIANEFSDRC TPSVISFGSK NRTIGVAAKN QQITHANNTV SNFKRFHGRA FNDPFIQKEK ENLSYDLVPL KNGGVGIKVM YMGEELFSV EQITAMLLTK LKETAENSLK KPVTDQVISV PSFFTDARR SVLDAAQIVG LNCLRLMNDM TAVALNYGIY KQDLPS

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Recombinant human HSPH1 protein SDS-PAGE gel 3µg reducing conditions Coomassie stain