

## Human Hexokinase 4 protein (Recombinant) (His-Tag) (STJP019670)

### GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

### PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol
Molecular Weight	54.3 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">2645</a>
Gene Symbol	<a href="#">GCK</a>
UniProt ID	<a href="#">HXK4_HUMAN</a>
Immunogen Region	1-465aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MLDDPAREMA AKKEKVEQIL AEFQLQEEDL KVMRRMQKE MDRGLRLETH EEASVKMLPT YVRSTPEGSE VGDFLSLDLG GTNFRVMLVK VGEGEEGQWS VKTKHQMYSI PEDAMTGTAELMLFDYISECI SDFLDKHQMK HKKLPLGFTF SFPVRHEDID KGILLNWTGK FKASGAEGNN VGLLRDAIK RRGDFEMDVV AMVNDTVA

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



15% SDS-PAGE (3µg)

SDS-PAGE gel of recombinant human Hexokinase 4 protein under reducing conditions  
Coomassie stained