

Human Hexokinase 2 protein (Recombinant) (His-Tag) (STJP019671)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Enzyme Activity
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol
Molecular Weight	104.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3099
Gene Symbol	HK2
UniProt ID	H XK2_HUMAN
Immunogen Region	1-917aa
Immunogen Sequence	MIASHLLAYF FTELNHDQVQ KVDQYLYHMR LSETLLEIS KRFRKEMEGK LGATTHPTAA VKMLPTFVRS TPDGTEHGEF LALDLGGTNF RVLWVKVTDN GLQKVMENQ IYAIPEIDMR GSGTQLFDHI AECLANFMDK LQIKDKKLPL GFTFSFPCHQ TKLDESFLVS WTKGFKSSGV EGRDVVALIR KAIQRRGDFD IDIVAV

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human Hexokinase 2 protein STJP019671 under reducing conditions with Coomassie blue stain