

Human IL-16 protein (Recombinant) (His-Tag) (STJP019686)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) containing 10% Glycerol
Molecular Weight	15.5 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3603
Gene Symbol	IL16
UniProt ID	IL16_HUMAN
Immunogen Region	502-631aa
Immunogen Sequence	MPDLNSSTD S AASASAASDV SVETAEATV CTVTLEKMSA GLGFSLEGGK GSLHGDKPLT INRIFKAAS EQSETVQPGD EILQLGGTAM QGLTRFEAWN IIKALPDGPV TIVIRRKSLQ SKETTAAGDS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE of recombinant human IL-16 protein under reducing conditions Coomassie blue stain