

Human IRF1 protein (Recombinant) (His-Tag) (STJP019689)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol
Molecular Weight	15 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3659
Gene Symbol	IRF1
UniProt ID	IRF1_HUMAN
Immunogen Region	1-114aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MPITRMRRP WLEMQINSNQ IPGLIWINKE EMIFQIPWKH AAKHGWDINK DACLFERSWAI HTGRYKAGEK EPDPKTWKAN FRCAMNSLPD IEEVKDQSRN KGSSAVRVYR MLPP

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE of recombinant human IRF1 protein under reducing conditions stained with Coomassie blue