

E.coli NusA protein (Recombinant) (No-Tag) (STJP019713)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4)
Molecular Weight	54 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Immunogen Region	1-495aa
Immunogen Sequence	MNKEILAVVE AVSNEKALPR EKIFEALESA LATATKKKYE QEIDVRVQID RKSGDFDTFR RWLVVDEVQ PTKEITLEAA RYEDESLNLG DYVEDQIESV TFDRIITTQTA KQVIVQKVR E AERAMVVDQF REHEGEIITG VVKKVNDRNI SLDLGNNAEA VILREDMLPR ENFRPGDRVR GVLVSVRPEA RGAQLFVTRS KPEMLIELFR IEVPEIGEEV IEIKAAAR

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant E.coli NusA protein 3µg under reducing conditions Coomassie blue stain