

Human PGP9.5/UCHL1 protein (Recombinant) (No-Tag) (STJP019723)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Enzyme Activity
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 2mM EDTA
Molecular Weight	24.8 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7345
Gene Symbol	UCHL1
UniProt ID	UCHL1_HUMAN
Immunogen Region	1-223aa
Immunogen Sequence	MLQKPMEINP EMLNKVLSRL GVAGQWRFVD VLGLEEEESLG SVPAPACALL LLFPLTAQHE NFRKKQIEEL KGQEVSPKVY FMKQTIGNSC GTIGLIHAVA NNQDKLGFED GSVLKQFLSE TEKMSPEDRA KCFEKNEAIQ AAHDAVAQEG QCRVDDKYNF HFILFNNVDG HLYELDGRMP FPNVHGASSE DTLLKDAAKV CREFTEREQG EVRFSAVALC KAA

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human PGP9.5 UCHL1 protein 3µg reducing conditions Coomassie blue stain