

NKMAX Marker (8 bands, 15-150 kDa) protein (Recombinant) (No-Tag) (STJP019726)

GENERAL INFORMATION


Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Liquid in 12mM Tris-HCl buffer (pH 6.8) 5% Glycerol, 0.4% SDS 2.88mM 2-mercaptoethanol, 0.02% bromophenol blue
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



15% SDS PAGE (5 μl)

NKMAX Protein Marker 8 bands SDS-PAGE gel under reducing conditions stained with Coomassie blue