

Human IPP-POZ (1-157aa) protein (Recombinant) (No-Tag) (STJP019728)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 10mM HEPES buffer (pH 7.4) containing 25mM NaCl
Molecular Weight	17.3 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3652
Gene Symbol	IPP
UniProt ID	IPP_HUMAN
Immunogen Region	1-157aa
Immunogen Sequence	MANEDCPKAA DSPFSSDKHA QLILAQINKM RNGQHFCDVQ LQVGQESFKA HRLVLAASSP YFAALFTGGM KESSKDVVPI LGIEAGIFQI LLDFIYTGIV NIGVNNVQEL IIAADMLQLT EVVHLCCEFL KGQIDPLNCI GIFQFSEQIA CHDLLEF

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human IPP-POZ 1-157aa protein under reducing conditions stained with Coomassie blue