

Human Cyclophilin F/PPIF protein (Recombinant) (His-Tag) (STJP019731)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Enzyme Activity
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 7.5) containing 1mM DTT, 10% Glycerol
Molecular Weight	21.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10105
Gene Symbol	PPIF
UniProt ID	PPIF_HUMAN
Immunogen Region	30-207aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MCSKGSQDPS SSSSSGNPLV YLDVDANGKP LGRVVLELKA DVVPKTAENF RALCTGEKGF GYKGFSTFHRV IPSFMCQAGD FTNHNGTGGK SIYGSRFPDE NFTLKHVGPV VLSMANAGPN TNGSQFFICT IKTDWLDGKH VVFGHVKEGM DVVKKIESFG SKSGRTSKKI VITDCGQLS

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing recombinant human Cyclophilin F PPIF protein under reducing conditions