

Human Cyclophilin H/PPIH protein (Recombinant) (No-Tag) (STJP019733)

GENERAL INFORMATION

| | |
|---------------|--------------------------|
| Product Type | Proteins |
| Applications | SDS-PAGE/Enzyme Activity |
| Host / Source | E.coli |

PRODUCT PROPERTIES


| | |
|---------------------|--|
| Concentration | 1 mg/mL |
| Formulation | Liquid in phosphate-Buffered Saline (pH 7.4) containing 10% Glycerol |
| Molecular Weight | 19.2 kDa |
| Storage Instruction | For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| | |
|--------------------|--|
| Gene ID | 10465 |
| Gene Symbol | PPIH |
| UniProt ID | PPIH_HUMAN |
| Immunogen Region | 1-177aa |
| Immunogen Sequence | MAVANSSPVN PVVFFDVSIG GQEVGRMKIE LFADVVPKTA ENFRQFCTGE FRKDGVPIGY KGSTFHRVIK DFMIQGGDFV NGDGTGVASI YRGPFADEFN KLRHSAPGLL SMANSGPSTN GCQFFITCSK CDWLDGKHVV FGKIIDGLLV MRKIENVPTG PNNKPKLPVV ISQCGEM |

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human Cyclophilin H PPIH protein under reducing conditions