

## Human PRMT1 protein (Recombinant) (His-MBP-Tag) (STJP019738)

### GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Enzyme Activity
Host / Source	E.coli

### PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 40mM Tris-HCl buffer (pH 8.0) containing 100mM NaCl, 4mM MgCl <sub>2</sub> , 2mM DTT, 40% Glycerol
Molecular Weight	84.2 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">3276</a>
Gene Symbol	<a href="#">PRMT1</a>
UniProt ID	<a href="#">ANM1_HUMAN</a>
Immunogen Region	1-353aa
Immunogen Sequence	

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human PRMT1 protein 3µg reducing conditions Coomassie blue stain