

Human RAB5A protein (Recombinant) (No-Tag) (STJP019745)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM MES buffer (pH 6.0) containing 10% Glycerol, 0.1mM PMSF
Molecular Weight	23 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5868
Gene Symbol	RAB5A
UniProt ID	RAB5A_HUMAN
Immunogen Region	1-215aa
Immunogen Sequence	MASRGATRPN GPNTGNKICQ FKLVLGESA VGKSSLVLRV VKGQFHEFQE STIGAAFLTQ TVCLDDTTVK FEIWDTAGQE RYHSLAPMY RGAQAAIVVY DITNEESFAR AKNWVKELQR QASPNIVIAL SGNKADLANK RAVDFQEAQS YADDNSLLFM ETSAKTSMNV NEIFMAIAKK LPKNPQNPG ANSARGRGVD LPTPTQPTRN QCCSN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human RAB5A protein under reducing conditions