

Human RAR alpha/NR1B1 protein (Recombinant) (His-Tag) (STJP019748)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid In Phosphate-Buffered Saline (pH 7.4) containing 10% Glycerol, 1mM DTT, 0.1M NaCl
Molecular Weight	14 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5914
Gene Symbol	RARA
UniProt ID	RARA_HUMAN
Immunogen Region	68-173aa
Immunogen Sequence	MSEEIVSPSP SPPPLPRIYK PCFVCQDKSS GYHYGVSACE GCKGFFRRSI QKNMVYTCHR DKNCIINKVT RNRCQYCRLQ KCFEVGMSKE SVRNDNRNKKK KEVPKPE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human RAR alpha NR1B1 protein under reducing conditions stained with Coomassie blue