

Human Retinol Binding 4/RBP4 protein (Recombinant) (No-Tag) (STJP019749)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) containing 10% Glycerol
Molecular Weight	21 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5950
Gene Symbol	RBP4
UniProt ID	RET4_HUMAN
Immunogen Region	19-201aa
Immunogen Sequence	MERDCRVSSF RVKENFDKAR FSGTWYAMAK KDPEGLFLQD NIVAEFSVDE TGQMSATAKG RVRLNNWDV CADMVGTFD TEDPAKFKMK YWGVASFQK GNDDHWIVDT DYDTYAVQYS CRLLNLDGTC ADSYFVFSR DPNGLPPEAQ KIVRQRQEEL CLARQYRLIV HNGYCDGRSE RNLL

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing recombinant human RBP4 protein 3µg under reducing conditions
Coomassie blue stain