


Human SETD7 protein (Recombinant) (No-Tag) (STJP019756)

GENERAL INFORMATION	
Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli
PRODUCT PROPERTIES	
Concentration	1 mg/mL
Formulation	Liquid in 50mM Tris-HCl buffer (pH 7.5) containing 0.2M NaCl, 5mM DTT, 20%Glycerol
Molecular Weight	40.7 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	80854
Gene Symbol	SETD7
UniProt ID	SETD7_HUMAN
Immunogen Region	1-366aa
Immunogen Sequence	MDSDDDEMVEE AVEGHLDDDG LPHGFCTVTY SSTDRFEGNF VHGEKNRGK FFFFDGSTLE GYYVDDALQG QGVYTYEDGG VLQGTYYDGE LNGPAQEYDT DGRILIFKGQY KDNIRHGVCW IYYPDGGSLV GEVNEDGEMT GEKIAYVYPD ERTALYGKFI DGEMIEGKLA TLMSTEEGRP HFELMPGNSV YHFDKSTSSC ISTNALLPDP YESERVYVAE SLISSAGE
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human SETD7 protein 3µg under reducing conditions with Coomassie blue stain