

Human Sirtulin 6/SIRT6 protein (Recombinant) (His-Tag) (STJP019759)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 7.5) containing 20% Glycerol, 1mM DTT, 0.1M NaCl, 0.1mM PMSF
Molecular Weight	41 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51548
Gene Symbol	SIRT6
UniProt ID	SIR6_HUMAN
Immunogen Region	1-355aa
Immunogen Sequence	MSVNYAAGLS PYADKGGKGL PEIFDPPEEL ERKVWELARL VWQSSSVVFH TGAGISTASG IPDFRGPBGV WTMEERGLAP KFDTTFESAR PTQTHMALVQ LERVGLLRFL VSNQVDGLHV RSGFPRDKLA ELHGNMFVEE CAKCKTQYVR DTVVGTMLGK ATGRLCTVAK ARGLRACRGE LRDITLDWED SLPDRDLALA DEASRN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human SIRT6 protein under reducing conditions stained with Coomassie blue