

Human SMAC/Diablo protein (Recombinant) (T7-Tag) (STJP019761)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 7.5)
Molecular Weight	22 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	56616
Gene Symbol	DIABLO
UniProt ID	DBLOH_HUMAN
Immunogen Region	56-239aa
Immunogen Sequence	MASMTGGQQM GRGSMAPVIA QKSEPHSLSS EALMRRVAVSL VTDSTSTFLS QTTYALIEAI TEYTKAVYTL TSLYRQYTSI LGKMNSEED EVWQVIIGAR AEMTSKHQY LKLETTWMTA VGLSEMAAEA AYQTGADQAS ITARNHIQLV KLQVEEVHQL SRKAETKLA E AQIEELRQKT QEEGEERAES EQEAYLRD

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human SMAC/Diablo protein under reducing conditions
Coomassie blue stain