

Human SNAPIN protein (Recombinant) (His-Tag) (STJP019763)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 5mM DTT, 2mM EDTA, 0.2 M NaCl2, 40% Glycerol
Molecular Weight	17 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23557
Gene Symbol	SNAPIN
UniProt ID	SNAPN_HUMAN
Immunogen Region	1-136aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MAGAGSAAVS GAGTPVAGPT GRDLFAEGLL EFLRPAVQQL DSHVHAVRES QVELREQIDN LATELCRINE DQKVALDLDV YVKLLNARR RVVLVNNILQ NAQERLRRLLN HSVAKETARR RAMLDGIYP PGSPGK

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



15% SDS-PAGE (3µg)

SDS-PAGE gel of recombinant human SNAPIN protein 3µg under reducing conditions Coomassie blue stain