

Human alpha-Synuclein 112/NACP112 protein (Recombinant) (No-Tag) (STJP019769)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 7.5) containing 0.1M NaCl
Molecular Weight	11.3 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6622
Gene Symbol	SNCA
UniProt ID	SYUA_HUMAN
Immunogen Region	1-112aa
Immunogen Sequence	MDVFMKGLSK AKEGVVAAAE KTKQGVAAEA GKTKEGVLYV GSKTKEGVVH GVATVAEKT EQVTNVGGAV VTGVTAVAQK TVEGAGSIAA ATGFVKKDQL GKEGYQDYEP EA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human alpha-Synuclein 112 NACP112 protein under reducing conditions