

Human alpha-Synuclein (61-140aa) protein (Recombinant) (No-Tag) (STJP019770)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 7.5) containing 0.1M NaCl
Molecular Weight	8.4 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6622
Gene Symbol	SNCA
UniProt ID	SYUA_HUMAN
Immunogen Region	61-140aa
Immunogen Sequence	MEQVTNVGGAV VTGVTAVAQK TVEGAGSIAA ATGFVKKDKL GKNEEGAPQE GILEDMPVDP DNEAYEMPSE EGYQDYEP EA

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing recombinant human alpha-Synuclein 61-140aa protein 3µg reducing conditions Coomassie blue stain