

## Human alpha-Synuclein (96-140aa) protein (Recombinant) (No-Tag) (STJP019771)

### GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

### PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 7.5) containing 0.1M NaCl
Molecular Weight	5.2 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">6622</a>
Gene Symbol	<a href="#">SNCA</a>
UniProt ID	<a href="#">SYUA_HUMAN</a>
Immunogen Region	96-140aa
Immunogen Sequence	MKKDQL GKNEEGAPQE GILEDMPVDP DNEAYEMPSE EGYQDYEP EA

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human alpha-Synuclein 96-140aa under reducing conditions  
Coomassie blue stain