

Human SNAP23 protein (Recombinant) (No-Tag) (STJP019780)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0)
Molecular Weight	23.3 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8773
Gene Symbol	SNAP23
UniProt ID	SNP23_HUMAN
Immunogen Region	1-211aa
Immunogen Sequence	MDNLSSEIIQ QRAHQITDES LESTRILGL AIESQDAGIK TITMLDEQKE QLNRIEEGLD QINKDMRETE KTLTELNKCC GLCVPCNRT KNFESGKAYK TTWGDGGENS PCNVVSKQPG PVTNGQLQQP TTGAASGGYI KRITNDARED EMEENLTQVG SILGNLKDMA LNIGNEIDAQ NPQIKRITDK ADTNRDRIDI ANARAKKLID S

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel showing recombinant human SNAP23 protein under reducing conditions
Coomassie blue stain