

Human Ubc9/UBE2I protein (Recombinant) (No-Tag) (STJP019794)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 50mM HEPES buffer (pH 7.4) containing 150mM NaCl, 1mM DTT, 10%Glycerol
Molecular Weight	18 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7329
Gene Symbol	UBE2I
UniProt ID	UBC9_HUMAN
Immunogen Region	1-158aa
Immunogen Sequence	MSGIALSRLA QERKAWRKDH PFGFVAVPTK NPDGTMNLMN WECAIPGKKK TPWEGGLFKL RMLFKDDYPS SPPKCKFEPP LFHPNVYPSPG TVCLSILEED KDWRPAITIK QILLGIQELL NEPNIQDPAQ AEAYTIYCQN RVEYEKRVRA QAKKFAPS

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human Ubc9 UBE2I protein under reducing conditions