

Human Ubiquitin/RPS27A protein (Recombinant) (No-Tag) (STJP019795)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 50mM HEPES buffer (pH 7.4) containing 150mM NaCl, 1mM DTT, 10% Glycerol
Molecular Weight	8.6 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6233
Gene Symbol	RPS27A
UniProt ID	RS27A_HUMAN
Immunogen Region	1-76aa
Immunogen Sequence	MQIFVKLTG KTITLEVEPS DTIENVKAKI QDKEGIPPDQ QRLIFAGKQL EDGRTLSDYN IQKESTLHLV LRLRGG

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



15% SDS-PAGE (3µg)
 SDS-PAGE gel of recombinant human Ubiquitin RPS27A protein under reducing conditions stained with Coomassie blue