

Human Ubch7/UBE2L3 protein (Recombinant) (No-Tag) (STJP019796)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 50mM HEPES buffer (pH 7.4) containing 150mM NaCl, 1mM DTT, 10% Glycerol
Molecular Weight	17.9 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7332
Gene Symbol	UBE2L3
UniProt ID	UB2L3_HUMAN
Immunogen Region	1-154aa
Immunogen Sequence	MAASRRLMKE LEEIRKCGMK NFRNIQVDEA NLLTWQGLIV PDNPPYDKGA FRIEINFP AE YPFKPPKITF KTKIYHPNID EKGQVCLPVI SAENWK PATK TDQVIQSLIA LVNDPQPEHP LRADLAE EYS KDRKKFCKNA EEFTKKYGEK RPVD

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing recombinant human Ubch7 UBE2L3 protein 3µg under reducing conditions stained with Coomassie blue