

Human VHL protein (Recombinant) (His-Tag) (STJP019804)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in phosphate-buffered Saline (pH 7.4) containing 1mM DTT, 2mM EDTA
Molecular Weight	19.2 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7428
Gene Symbol	VHL
UniProt ID	VHL_HUMAN
Immunogen Region	1-154aa
Immunogen Sequence	MPRAENWDE AEVGAEAGV EEYGPEDDG EESGAEESGP EESGPEELGA EEEMEAGRPR PVLRSVNSRE PSQVIFCNRS PRVLPVWLN FDGEPQPYPT LPPGTGRRRH SYRGHLWLFV DAGTHDGLLV NQTELVFVPSL NVDGQPIFAN ITLP

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Recombinant human VHL protein Cat No. STJP019804 antibody reagent