

Human 14-3-3 zeta protein (Recombinant) (His-Tag) (STJP019812)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4)
Molecular Weight	32 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7534
Gene Symbol	YWHAZ
UniProt ID	1433Z_HUMAN
Immunogen Region	1-245aa
Immunogen Sequence	MDK NELVQKAKLA EQAERYDDMA ACMKSVTEQG AELSNEERNL LSVAYKNVVG ARSSWRVVS SIEQKTEGAE KKQQMAREYR EKIE TELRDI CNDVLSLLEK FLIPNASQAE SKVFYLMKMG DYYRYLAEVA AGDDKKGIVD QSQQAYQAEF EISKEMQPT HPIRLGLALN FSVFYEILN SPEKAC

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human 14-3-3 zeta protein under reducing conditions stained with Coomassie blue