

Human Siglec-4/MAG/GMA protein (Recombinant) (C-hFc&His) (STJP019813)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4099
Gene Symbol	MAG
UniProt ID	MAG_HUMAN
Immunogen Region	Gly20-Pro516
Immunogen Sequence	GHWGAWMPSSISAFEGTCVS IPCRFDFPDELRPVVHGVW YFNSPYPKNYPPVVKSRQTQ VVHESFQGRSRLLDLGLRN CTLLSNVSPELGGKYYFRG DLGGYNQYTFSEHSVLDIVN TPNIVPPEVVAGTEVEVSC MVPDNCPELRPELSWLGHEG LGEPAVLGRLREDEGTWVQV SLLHFVPTREANGHRLGCQA SFPNTTLQFEGYASMDVKYP PVIVEMNSSVEAIEGSHVS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---