

Human Cripto-1/TDGF1 protein (Recombinant) (C-His) (STJP019814)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.Contact us for customized product form or formulation.
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6997
Gene Symbol	CRIPTO
UniProt ID	TDGF1_HUMAN
Immunogen Region	Leu31-Thr172
Immunogen Sequence	LGHQEFARPSRGYLAFRDDS IWPQEEP AIRPRSSQRVPPM GIQHSKELNRTCCLNGGTCM LGSFCACPPSFYGRNCEHDV RKENCGSVPHDTWLPKKCSL CKCWHGQLRCFPQAFLPGCD GLVMDEHLVASRTELPSPA RT

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---