

Human B7-DC/PD-L2/CD273 protein (Recombinant) (C-His&Avi) {Biotin} (STJP019816)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	80380
Gene Symbol	PDCD1LG2
UniProt ID	PD1L2_HUMAN
Immunogen Region	Leu20-Thr220
Immunogen Sequence	LFTVTVPKELYIIHGSNVT LECNFDTGSHVNLGAIASL QKVENDTSPHRERATLLEEQ LPLGKASFHIPQVQVRDEGQ YQCIIIYGVAWDYKYLTLKV KASYRKINTHILKVPETDEV ELTCQATGYPLAEVSWPNVS VPANTSHSRTPEGLYQVTSV LRLKPPPGRNFSCVFWNTHV RELTASIDLQSQMEPRTHP T

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---