

Human B7-H5/Gi24/VISTA protein (Recombinant) (C-His) (STJP019819)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	64115
Gene Symbol	VSIR
UniProt ID	VISTA_HUMAN
Immunogen Region	Phe33-Ala194
Immunogen Sequence	FKVATPYSLYVCPGQNVTL TCRLGPDVKGHDTVYKWTW YRSSRGEVQTCSERRPIRNL TFQDLHLHHGGHQAANTSHD LAQRHGLESASDHGNSIT MRNLTLDSGLYCCLVVEIR HHHSEHRVHGAMELQVQTGK DAPSNCVVYPSSSQDSENIT AA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---