

Human SLAMF7/CRACC/CD319 protein (Recombinant) (C-mFc) (STJP019822)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	57823
Gene Symbol	SLAMF7
UniProt ID	SLAF7_HUMAN
Immunogen Region	Ser23-Met226
Immunogen Sequence	SGPVKELVGSVGGAVTFPLK SKVKQVDSIVWTFNTTPLVT IQPEGGTIIVTQNRNRERVD FPDGGYSLKLSKLKKNDSGI YYVGIYSSSLQQPSTQEYVL HVYEHLSPKPVMTMGLQSNKN GTCVTNLTCMEHGEEVDVIY TWKALGQAANESHNGSILPI SWRWGESDMTFICVARNPVS RNFSSPILARKLCEGAADDP DSSM

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---