

## Human SLAMF7/CRACC/CD319 protein (Recombinant) (C-hFc) (STJP019823)

| GENERAL INFORMATION    |                                                                                                                                                                                                                                   |
|------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type           | Proteins                                                                                                                                                                                                                          |
| Host / Source          | HEK293 cells                                                                                                                                                                                                                      |
| PRODUCT PROPERTIES     |                                                                                                                                                                                                                                   |
| Formulation            | Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.                                                                                                 |
| Storage Instruction    | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                                                                                                                                    |
| TARGET INFORMATION     |                                                                                                                                                                                                                                   |
| Gene ID                | <a href="#">57823</a>                                                                                                                                                                                                             |
| Gene Symbol            | <a href="#">SLAMF7</a>                                                                                                                                                                                                            |
| UniProt ID             | <a href="#">SLAF7_HUMAN</a>                                                                                                                                                                                                       |
| Immunogen Region       | Ser23-Met226                                                                                                                                                                                                                      |
| Immunogen Sequence     | SGPVKELVGSVGGAVTFPLK SKVKQVDSIVWTFNTTPLVT IQPEGGTIIVTQNRNRERVD<br>FPDGGYSLKLSKLKKNDSGI YYVGIYSSSLQQPSTQEYVL HVYEHLSPKPVMTMGLQSNKN<br>GTCVTNLTCCMEHGEEVDVIY TWKALGQAANESHNGSILPI SWRWGESDMTFICVARNPVS<br>RNFSSPILARKLCEGAADDP DSSM |
| ADDITIONAL INFORMATION |                                                                                                                                                                                                                                   |
| Note                   | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.                                                                                                           |