

Human SIRP-gamma/CD172g protein (Recombinant) (C-His&Avi) (STJP019824)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	55423
Gene Symbol	SIRPG
UniProt ID	SIRPG_HUMAN
Immunogen Region	Glu29-Pro360
Immunogen Sequence	EELQMIQPEKLLLVTVGKTA TLHCTVTSLLPVGPVLWFRG VGPGRELIYNQKEGHFPRVT TVSDLTKRNNMDFSIIRISSI TPADVGTYYCVKFRKGGSPEN VEFKSGPGTEMALGAKPSAP VVLGPAARTTPEHTVSFTCE SHGFSPRDITLKWFKNGNEL SDFQTNVDPTGQSVAYSIRS TARVVLDPWDVRSQVICEVA HVTLQGDPLRGTANLSEAIR VPPTLEVTQQPMRVGNQVN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---