

Human DLL4 protein (Recombinant) (C-His&Avi) {Biotin} (STJP019826)

GENERAL INFORMATION	
Product Type	Proteins
Host / Source	HEK293 cells
PRODUCT PROPERTIES	
Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS, 200 mM L-arginine (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	54567
Gene Symbol	DLL4
UniProt ID	DLL4_HUMAN
Immunogen Region	Ser27-Pro524
Immunogen Sequence	SGVFLQLQLQEFINERGLAS GRPCEPGCRTFFRVCLKHFQ AVVSPGPCTFGTVSTPVLGT NSFAVRDDSSGGGRNPLQLP FNFTWPGTFSLIIEAWHAPG DDLRPEALPPDALISKIAIQ GSLAVGQNWLLDEQTSTLTR LRYSYRVICSDNYGDNCSR LCKKRNDHFGHYVCQPDGNL SCLPGWTGEYCQQPICLSGC HEQNGYCSKPAECLCRPGWQ GRLCNECIPHNGCRHGTC
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.