

Human TREM-1/CD354 protein (Recombinant) (C-His) (STJP019827)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54210
Gene Symbol	TREM1
UniProt ID	TREM1_HUMAN
Immunogen Region	Ala21-Arg200
Immunogen Sequence	ATKLTEEKYELKEGQTLDVK CDYTLKFASSQKAWQIIRD GEMPKTLACTERPSKNSHPV QVGRRIILEDYHDHGLLRVRM VNLQVEDSGLYQCVIYQPPK EPHMLFDRIRLVVTGFGSGT PGSNENSTQNVIKIPPTTTK ALCPLYTSPRTVTQAPPKST ADVSTPDSEINLTNVTDIIR

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---