

Human TREM2 protein (Recombinant) (C-hFc) (STJP019828)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS, 100 mM Arginine (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54209
Gene Symbol	TREM2
UniProt ID	TREM2_HUMAN
Immunogen Region	His19-Ser174
Immunogen Sequence	HNTTVFQGVAGQSLQVSCP Y DSMKHWGRRKAWCRQLGEKG PCQRVVSTHNLWLLSFLRRW NGSTAITDDTLGGTLTITLR NLQPHDAGLYQCQSLHGSEA DTLRKVLVEVLADPLDHRDA GDLWFPGESESFEDAHVEHS ISRSLEGEIPFPPTS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---