

Human B7-H1/PD-L1/CD274 protein (Recombinant) (C-His) (STJP019830)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	29126
Gene Symbol	CD274
UniProt ID	PD1L1_HUMAN
Immunogen Region	Phe19-Ile226
Immunogen Sequence	FTVTVPKDLVVEYGSNMTI ECKFPVEKQLDLAALIVYWE MEDKNIIQFVHGEECLKVQH SSYRQRRARLLKDQLSLGNAA LQITDVKLQDAGVYRCMISY GGADYKRITVKVNAPYNKIN QRILVDPVTSEHELTCAE GYPKAEVIWTSSDHQVLSGK TTTTNSKREEKLFNVTSTLR INTTTNEIFYCTFRRLDPEE NHTAELVI

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---