

Human DNAM-1/CD226 protein (Recombinant) (C-hFc) (STJP019835)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10666
Gene Symbol	CD226
UniProt ID	CD226_HUMAN
Immunogen Region	Glu19-Asn247
Immunogen Sequence	EEVLWHTSVPFPAENMSLECV YPSMGILTQVEWFKIGTQQD SIAIFSPHGMVIRKPYAER VYFLNSTMASNNMTLFFRNA SEDDVGYYSCLYTPQGTW QKVIQVVQSDSFEAAVPSNS HIVSEPGKNVTLTCQPQMTW PVQAVRWEEKIQPRQIDLLTY CNLVHGRNFTSKFPRQIVSN CSHGRWSVIVIPDVTVSDSG LYRCYLQASAGENETFMRL TVAEGKTDN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---