

Human LILRB2/ILT-4/CD85d protein (Recombinant) (C-His&Avi) {Biotin} (STJP019842)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10288
Gene Symbol	LILRB2
UniProt ID	LIRB2_HUMAN
Immunogen Region	Gln22-His458
Immunogen Sequence	QTGTIPKPTLWAEPSVITQ GSPVTLSCQGSLEAQEYRLY REKKSASWITRIRPELVKNG QFHIPSITWEHTGRYGCQYY SRARWSELSDDLVLVMTGAY PKPTLSAQSPVVTSGGRVT LQCESQVAFGGFILCKEGED EHPQCLNSQPHARGSSRAIF SVGPVSPNRRWSHRCYGYDL NSPYVWSSPSDLELLVPGV SKKPSLSVQPGPVMAPGESL TLQCVSDVG YDRFVLYKEG

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---