

Human TACTILE/CD96 protein (Recombinant) (C-mFc) (STJP019846)

GENERAL INFORMATION	
Product Type	Proteins
Host / Source	HEK293 cells
PRODUCT PROPERTIES	
Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	10225
Gene Symbol	CD96
UniProt ID	TACT_HUMAN
Immunogen Region	Val22-Met503
Immunogen Sequence	VWEKTVNTEENVYATLGSDV NLTCQTQTVGFFVQMOWSKV TNKIDLIAYVHPQYGFYCA Y GRPCESLVTFETEPENGSKW TLHLRNMSCSVSGRYECMLV LYPEGIQTKIYNLLIQTHVT ADEWNSNHTIEIEINQTL E PCFQNSSSKISSEFTYAWSV EDNGTQETLISQNHLSNST LLKDRVKLGTDYRLHLSVPQ IFDDGRKF SCHIRVGP NKIL RSSTTVKVF AKPEIPVIVE
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.