

Human NKp46/NCR1/CD335 protein (Recombinant) (C-mFc) (STJP019848)

GENERAL INFORMATION	
Product Type	Proteins
Host / Source	HEK293 cells
PRODUCT PROPERTIES	
Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	9437
Gene Symbol	NCR1
UniProt ID	NCTR1_HUMAN
Immunogen Region	Gln22-Asn255
Immunogen Sequence	QQQTLPKPFIWAEPHFMPVK EKQVTICCGNYGAVEYQLH FEGSLFAVDRPKPPERINKV KFYIPDMNSRMAGQYSCIYR VGELWSEPSNLLDLVVTEMY DTPTLSVHPGPEVISGEKVT FYCRDLDATSMFLLKKEGRS SHVQRGYGVQAEFFPLGPVT TAHRGTYRCFGSYNNHAWSF PSEPVKLLVTGDIENTSLAP EDPTFPADTWGTYLLTTETG LQKDHALWDHTAQN
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.