

Human Siglec-5/CD170 protein (Recombinant) (C-His&Avi) {Biotin} (STJP019850)

GENERAL INFORMATION	
Product Type	Proteins
Host / Source	HEK293 cells
PRODUCT PROPERTIES	
Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS, 350 mM NaCl (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	8778
Gene Symbol	SIGLEC5
UniProt ID	SIGL5_HUMAN
Immunogen Region	Glu17-Thr434
Immunogen Sequence	EKPVYELQVQKSVTVQEGLC VLVPCSFSPWRSWYSSPPL YVYWFRDGEIPYYAEVVATN NPDRRVKPETQGRFRLLGDV QKKNCSLSIGDARMEDTGSY FFRVERGRDVKYSYQQNKLN LEVTALIEKPDHIFLEPLES GRPTRLSCSLPGSCEAGPPL TFSWTGNALSPLDPETTRSS ELTLTPRPEDHGTNLTCQMK RQGAQVTTERTVQLNVSYAP QTITIFRNGIALEILQNTS
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.