

Human Siglec-5/CD170 (Primary Amine Labeling) protein (Recombinant) (C-hFc) {Biotin} (STJP019851)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8778
Gene Symbol	SIGLECS
UniProt ID	SIGL5_HUMAN
Immunogen Region	Glu17-Thr434
Immunogen Sequence	EKPVYELQVQKSVTVQEGLC VLVPCSFSPWRSWYSSPPL YVYWFRDGEIPYYAEVVATN NPDRRVKPETQGRFRLLGDV QKKNCSLSIGDARMEDTGSY FFRVERGRDVKYSYQQNKLN LEVTALIEKPDHIFLEPLES GRPTRLSCLPGSCEAGPPL TFSWTGNALSPLDPETTRSS ELTLTPRPEDHGTNLTCQMK RQGAQVTTERTVQLNVSYAP QTITIFRNGIALEILQNTS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---