

Human CD24 (Primary Amine Labeling) protein (Recombinant) (C-mFc) {Biotin} (STJP019876)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	100133941
Gene Symbol	CD24
UniProt ID	CD24_HUMAN
Immunogen Region	Ser27-Gly59
Immunogen Sequence	SETTTGTSSNSSQSTSNSGL APNPTNATTKAAG

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---