

Human 5 NA-Nucleotidase/NT5E/CD73 protein (Recombinant) (C-hFc) (STJP019883)

GENERAL INFORMATION	
Product Type	Proteins
Host / Source	HEK293 cells
PRODUCT PROPERTIES	
Formulation	Lyophilised from 0.22 M filtered solution in 20 mM Tris, 120 mM NaCl (pH 7.5). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	4907
Gene Symbol	NT5E
UniProt ID	5NTD_HUMAN
Immunogen Region	Trp27-Lys547
Immunogen Sequence	WELTILHTNDVHSRLEQTSE DSSKCVNASRCMGGVARLFT KVQQIRRAEPNVLLLDAGDQ YQGTIWFTVYKGA EVAHFMN ALRYDAMALGNHEFDNGVEG LIEPLLKEAKFPILSANIKA KGPLASQISGLYPYKVL PV GDEVVGIVGYTSKETPFSLN PGTNLVFEDEITALQPEVDK LKT LNVNKI IALGHSGFEMD KLIAQKVRGVDVVVGHSNT FLYTGNPPSKEVPAGKYPF
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.