

## Human CEACAM6/CD66c protein (Recombinant) (C-His&Avi) (STJP019884)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

### PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">4680</a>
Gene Symbol	<a href="#">CEACAM6</a>
UniProt ID	<a href="#">CEAM6_HUMAN</a>
Immunogen Region	Lys35-Gly320
Immunogen Sequence	KLTIESTPFNVAEGKEVLLL AHNLPQNRIGYSWYKGERVD GNSLIVGYVIGTQQATPGPA YSGRETIYPNASLLIQNVTQ NDTGFYTLQVIKSDLVNEEA TGQFHVYPELPKPSISSNNS NPVEDKDAVAFTCEPEVQNT TYLWWWVNGQSLPVSRLQLS NGNMTLTLVSKRNDAGSYE CEIQNPASANRSDPVTNLVL YGPDVPTISPSKANYRPGEN LNLSCHAASNPPAQYSWFI

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---