

Human CEACAM6/CD66c protein (Recombinant) (C-His&Avi) {Biotin} (STJP019885)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4680
Gene Symbol	CEACAM6
UniProt ID	CEAM6_HUMAN
Immunogen Region	Lys35-Gly320
Immunogen Sequence	KLTIESTPFNVAEGKEVLLL AHNLPQNRIGYSWYKGERVD GNSLIVGYVIGTQQATPGPA YSGRETIYPNASLLIQNVQT NDTGFYTLQVIKSDLVNEEA TGQFHVYPELPKPSISSNNS NPVEDKDAVAFTCEPEVQNT TYLWWVNGQSLPVSPRLQLS NGNMTLTLLSVKRNDAGSYE CEIQNPASANRSDPVTNLVL YGPDVPTISPSKANYRPGEN LNLSCHAAASNPPAQYSWFI

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---