

Human CD3E&CD3G (Primary Amine Labeling) protein (Recombinant) (C-His) {Biotin} (STJP019887)

GENERAL INFORMATION	
Product Type	Proteins
Host / Source	HEK293 cells
PRODUCT PROPERTIES	
Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Immunogen Region	Asp23-Glu120
Immunogen Sequence	DGNEEMGGITQTPYKVSISG TTVILTCPQYPGSEILWQHN DKNIGGDEDDKNIGSDEDHL SLKEFSELEQSGYYVCYPRG SKPEDANFYLYLRARVCEQS IKGNHLVKVYDYQEDGSVLL TCDAEAKNITWFKDGKMIGF LTEDKKKWNLGSNAKDPRGM YQCKGSQNKSKPLQVYYRMC QN
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.